


CURAPROX CPS Chairside Box prime

25x CPS 06-011 pre-assembled on UHS 409 holder + 125 cap, 24x IAP

Codes

CURADEN Article No.	73345033
Customs Pos. No.	9603.2100
Country of Origin	Germany



	Content
Content Box	25x CPS 06 pre-assembled on UHS 409 holder incl. protective cap 25x CPS 07 pre-assembled on UHS 409 holder incl. protective cap 25x CPS 08 pre-assembled on UHS 409 holder incl. protective cap 25x CPS 09 pre-assembled on UHS 409 holder incl. protective cap 25x CPS 011 pre-assembled on UHS 409 holder incl. protective cap 24 x IAP Prime

Product profile interdental brushes CPS prime 06

	Overall length	Diameter Ø
CPS 06	19 mm	2.2 mm

	Material	Diameter Ø	Colour
Brush	Nylon	0.05 mm	transparent

	Material	Length	Diameter Ø	Colour
Shaft	Polypropylene	8.6 mm	2.3 mm	turquoise

	Material	Diameter Ø
Core	Stainless steel	0.18 mm

CPS prime 07

	Overall length	Diameter Ø
CPS 07	20 mm	2.5 mm

	Material	Diameter Ø	Colour
Brush	Nylon	0.06 mm	transparent

	Material	Length	Diameter Ø	Colour
Shaft	Polypropylene	8.6 mm	2.3 mm	red

	Material	Diameter Ø
Core	Stainless steel	0.20 mm

SPECIFICATION SHEET

CURAPROX

CPS prime 08

	Overall length	Diameter Ø		
CPS 08	21.5 mm	3.2 mm		
	Material	Diameter Ø		Colour
Brush	Nylon	0.08 mm		transparent
	Material	Length	Diameter Ø	Colour
Shaft	Polypropylene	8.6 mm	2.3 mm	pink
	Material	Diameter Ø		
Core	Stainless steel	0.23 mm		

CPS prime 09

	Overall length	Diameter Ø		
CPS 09	21.5 mm	4.0 mm		
	Material	Diameter Ø		Colour
Brush	Nylon	0.08 mm		transparent
	Material	Length	Diameter Ø	Colour
Shaft	Polypropylene	8.6 mm	2.3 mm	yellow
	Material	Diameter Ø		
Core	Stainless steel	0.23 mm		

CPS prime 011

	Overall length	Diameter Ø		
CPS 011	22.0 mm	5.0 mm		
	Material	Diameter Ø		Colour
Brush	Nylon	0.08 mm		transparent
	Material	Length	Diameter Ø	Colour
Shaft	Polypropylene	8.6 mm	2.3 mm	lime green
	Material	Diameter Ø		
Core	Stainless steel	0.25 mm		

Product profil holder

	Material	weight	color
Holder	Tritan	2.5 g	turquoise, red, pink, yellow, lime green
Protective cap	Styrolbutadiencopolymer (SBC)	0.5 g	transparent

SPECIFICATION SHEET

CURAPROX

Product profile IAP

Material	length	diameter Ø	colours	weight
Polysulfon (PSU)	22,5mm	0.6 mm – 2.3 mm	turquoise, red, pink, yellow, lime wood	0.03 g

Product profil box

	Material	colour	weight	Dimensions (W x L x D)
Lid	Polystyrene (PS)	Transparent	171g	298x167x16mm
Inlay	Polystyrene (PS)	white	78g	288x157x38mm
box	Styrol / Butadien (SB)	blue	149g	288x157x38mm

Packaging

	Material	Dimensions (W x L x D)
Slipcase for box	Cardboard	344x173x46 mm

Units

	Content	Weight	Outer Measures (Lx W)
Box	1	642g	298x157x38mm
Consumer pack	1	814g	344x173x46 mm
Carton	20	18.6kg	58.5cm x 38.5cm x 39.5cm