

CURAPROX CPS 150 prime mixed set
with 1 x CPS 06, 07, 08, 09, 011 prime + UHS 409 holders

Codes

CURADEN Article No.	73340112
Customs Pos. No.	9603.2100
Country of Origin	Germany



Product profile brush

CPS prime 06

	Material	Length	Diameter Ø	Layer	Colour
Brush	Nylon	10mm	2.2mm	single	transparent
	Material	Length	Diameter Ø	Colour	
Shaft	Polypropylene	8.6mm	2.3mm	turquoise	
	Material	Length	Diameter Ø		
Core	Stainless steel	10mm	0.18mm		

CPS prime 07

	Material	Length	Diameter Ø	Layer	Colour
Brush	Nylon	10mm	2.5mm	single	transparent
	Material	Length	Diameter Ø	Colour	
Shaft	Polypropylene	8.6mm	2.3mm	red	
	Material	Length	Diameter Ø		
Core	Stainless steel	10mm	0.18mm		

CPS prime 08

	Material	Length	Diameter Ø	Layer	Colour
Brush	Nylon	12.5 mm	3.2mm	single	transparent
	Material	Length	Diameter Ø	Colour	
Shaft	Polypropylene	8.6mm	2.3mm	pink	
	Material	Length	Diameter Ø		
Core	Stainless steel	11.5mm	0.23 mm		

SPECIFICATION SHEET

CURAPROX

CPS prime 09

	Material	Length	Diameter Ø	Layer	Colour
Brush	Nylon	12.5 mm	4.0 mm	single	transparent
	Material	Length	Diameter Ø	Colour	
Shaft	Polypropylene	8.6 mm	2.3 mm	yellow	
	Material	Length	Diameter Ø		
Core	Stainless steel	12.0 mm	0.23 mm		

CPS prime 011

	Material	Length	Diameter Ø	Layer	Colour
Brush	Nylon	13 mm	5.0 mm	single	transparent
	Material	Length	Diameter Ø	Colour	
Shaft	Polypropylene	8.6 mm	2.3 mm	Light green	
	Material	Length	Diameter Ø		
Core	Stainless steel	12.0 mm	0.25 mm		

Product profile handle

Product profil hodler

	Material	weight	color
holder	Tritan	2.5 g	Turquoise, red, pink, yellow, light green
Protective cap	SBC	0.5 g	transparent

Consumer pack packaging

	Material	Weight	Dimensions (L x W)
Minigrip	Polyethylen (LDPE)	0.9 g	L 92mm x W 65mm
Enclosures	Paper	2.11 g	L 88mm x W 56mm

Units

	Content	Weight	Outer Measures (LxWxH)	Volume dm ³
Consumer pack include 4 brushes with holders	1	22.5 g	92mm x 65mm x 5mm	0.03
Carton	300	7.5 kg	39.5 x 30 x 22cm	18.06