

KY-826

*Professional trust,
Selected at a great value*

Ergonomically Designed, Comprehensive Care and Support




KUANG YEU MEDICAL Co., Ltd.
DESIGN & MANUFACTURING IN TAIWAN
TEL:+886 4 2350 8866
FAX:+886 4 2350 0099
www.kuang-yeu.com.tw
marketing2@kuang-yeu.com.tw
No.10, Industrial Zone 7, Xitun District, Taichung City



KY-826 (A)

- **Smart adjustment:** Supports preset treatment positions, making it easier for doctors to operate and for patients to get up.
- **Silent Drive System:** Eliminate mechanical noise, allowing you to focus on every detail of your treatment.
- **You can choose different operation modes** to suit your usage habits and improve work efficiency.



KY-826 (B)

- **The suspended pull-up cable system** delivers the most precise, non-invasive upright exercise solution. It prevents soreness from poor posture while distributing the weight evenly across your grip, allowing effortless and stress-free use.
- **The patients upper and lower body positions** do not touch the tubing, and there is no risk of instruments falling.
- **The rotating tray design** makes it easier to retrieve items, boosting your work efficiency.



KY-826 (C)

- Choose a mobile cart to expand your treatment range. The instrument surface remains fixed regardless of chair height adjustments, and its mobility allows you to maintain optimal working angles for various treatment tasks, significantly boosting your efficiency.



- Three-axis steering maintains design integrity with enhanced operational precision.
- Single-color LED with 4500K color temperature delivers optimal brightness and softness for extended use, protecting your eyes. Features a touch-sensitive, infinitely adjustable brightness range (8000-40000 LEX) to meet your lighting needs, with both manual and touch modes available.
- The chair features a convenient interlocking function that automatically turns off the light when you sit up and turns it on when you lie down.



- An intuitive functional panel design that enhances precise operation and reduces error risks;
- Features multi-angle phone holders for convenient operation and personalized usage.



- The funnel-shaped deep sputum basin is designed to prevent splashing when rinsing the mouth.
- The advanced infrared automatic sensing water supply system is user-friendly;
- The design of the sputum basin flushing and water supply sequence ensures the basin remains clean during use.

Best maneuverability, Presenting professional credibility



• From precision control to human-centered care, from efficient execution to safety protection, Guangyu Dental Chairs redefine oral healthcare standards through innovative technology, empowering dentists to focus on their expertise while maintaining patients unwavering trust.



• 90° pivot armrest design with snap-fit mechanism ensures unobstructed access to upper and lower examination positions.
• The optimal support ultra-short platform design provides maximum therapeutic space for your feet.




• Multifunctional foot pedal controller: Enables control of surgical lamp activation, dental chair elevation/tilt, mouth rinse, and chair repositioning, ensuring seamless contactless operation.



• Meeting the needs of different body types, providing optimal treatment positions

Customize to your needs Select accessories

• Multiple dental instrument applications, tailored to your clinical needs, delivering a customized professional operating system.



intraoral camera

Handmade, hand-sewn leather cushion for dental chairs

The fusion of comfort and durability in art

- Precision fit: Customized with individualized stitching to perfectly align with dental chair contours, ensuring maximum patient comfort.
- Exceptional durability: High-density sponge combined with imported medical-grade leather, offering tear resistance and color retention, with a lifespan tripled
- Premium detail craftsmanship: Hand-embroidered edges with precision, reinforced double-stitch construction, and warp-resistant design to prevent thread deformation
- Medical-grade PU/PVC leather: Easy to disinfect and resistant to bacterial growth
- Thickened high-rebound foam: Distributes pressure evenly, preventing collapse during prolonged sitting.

